GEOCHEMISTRY ANALYSIS: Geochemistry Equipment

Geochemistry Equipment Catalogue Carbon Analysis Equipment **C230 CARBON DETERMINATION**

The LECO C230 Carbon Determinator is a microprocessor based, software driven instrument for measurement of carbon and sulphur content in metals, ores, ceramics and other inorganic materials. The C230 uses an induction furnace and measures carbon by infrared absorption.

The infrared (IR) source consists of nichrome wire that is resistance-heated to 850°C. The IR source radiates visible energy as well as all wavelengths in the infrared spectrum. The description may use carbon dioxide; the same principle applies to sulfur detection.

Carbon dioxide absorbs IR energy at a precise wavelength within the IR spectrum. Energy from IR source is absorbed as the gas passes through the cell, preventing it from reaching the IR detector.

One IR cell is used both as a reference and for measurement. The total carbon, as carbon dioxide, is detected on a continuous and simultaneous basis.

The starting reference level, or baseline, for the detector is established by running oxygen through the cell to drive off residual atmospheric gases.

Specifications

Range (at 1 gram) C230C Carbon : 4 ppm to 3.5% C230CLC Carbon : 2 ppm to 0.5% **C230CHC** Carbon : 60 ppm to 6.0% Note: Reducing sample weight may extend n Precision C230C Carbon : 2 ppm or 0.5% RSI C230CLC Carbon : 1 ppm or 0.5% RSD* **C230CHC**

Carbon : 30 ppm or 0.5% RSD*



LECO C230 CARBON ANALYZER

	Result Readability		
	With PC : 15 digits of precision		
	Without PC : 0.1 ppm		
	Calibration		
	With PC : Multi-Point Linear		
	Without PC : Single-Point Linear		
range.	Analysis Time: 45 seconds (Nominal)		
	Sample size : 1 gram (Nominal)		
	Detection method : Sulphur State, Infrared Absorption,		
D*	Sulphur as Sulphur Dioxide		
	Chemical Reagents : Anhydrone Magnesium		

Perchlorate, Sodium Hydroxide on an inert base, Lecosorb, Platinized Silica. Cellulose

GEOCHEMISTRY ANALYSIS: Geochemistry Equipment

Geochemistry Equipment Catalogue Carbon Analysis Equipment

Specifications

Gas Required		Dimensions
Carrier Gas	: Oxygen 99.5% pure, 40 psi	Determinator
Pneumatic Gas	: Compressed Air, Nitrogen, or	Height :
	Argon, 40 psi (2.76 bars),	Width :
	Source must be oil & water free	
		Depth :
	on nominal analysis time)	Weight :
Measure	: 3 litres per analysis	Computer (PC b
Pneumatic	: 1 litre per analysis	Height :
D		0
Regulators		Width :
Oxygen	: 501-291 O ₂ Pressure Regulator	Depth :
Compressed Air	- 0	Weight:
		M 4 (DC1
Inert Gas	: 764-216 Inert Gas Regulator	Monitor (PC bas
Furnace	: Induction, 18 MHz, 2.2 kW	Height :
Furnace	• Induction, 10 MIL2, 2.2 KW	Width :
Data Transmit	: Included	Depth :
		Weight :
Data Storage		
With PC	: Limited only by hard drive	Electrical Power
	space	
Without PC	: 10 Samples Weight,	Determinator :
	50 Answers, 5 Methods	
Kev Pad	e o misirero, e methous	Operating Curr
With PC	: External	Stand-by Curren
With PC	: Internal Membrane /	Computer (PC b
without PC		• •
	External Optional	Monitor (PC bas
Display		
With PC	: 15.00 inch SVGA Monitor	Balance: 120V, =
Without PC	: Liquid Crystal (16x26 character)	Datalice, 120V,

ninator				
Height : 30.50 inches (77.5 cm)				
Width : 27.50 inches (70.0 cm)				
Depth : 23.50 inches (59.7 cm)				
Weight: 300 pounds				
uter (PC based system)				
Height : 17.00 inches (43.2 cm)				

8.00 inches (20.3 cm) 17.00 inches (43.2 cm) 29 pounds

sed system)

15.00 inches (38.1 cm) 14.00 inches (35.6 cm) 16.00 inches (40.6 cm)

35 pounds

Requirements

	Determinator : 230 \vec{V} , $\pm 10\%$, 50/60 Hz,
	Single phase, 15 amps max
	Operating Current: 12 Amps
	Stand-by Current : 4 Amps
	Computer (PC based): 115/230 V, ± 10%, 50/60
	Hz, 5/3 amps max
	Monitor (PC based): 90 264 V, 50/60 Hz,
or	1.6 amps max
naracter)	Balance: 120V, ± 10%, 50/60 Hz, 0.5 amps max

www.oiltestgroup.com